

# Recent developments in automated perimetry in glaucoma diagnosis and management

Chris A. Johnson, PhD

Recently, there have been several new developments in automated perimetry that have contributed to enhanced diagnosis and management of glaucoma. This paper will briefly review four of the latest advances in automated perimetry: (1) efficient test strategies that reduce the testing time and variability of automated perimetric testing, in particular, the Swedish Interactive Threshold Algorithm (SITA) and Tendency Oriented Perimetry (TOP) test strategies will be described; (2) Frequency Doubling Technology (FDT) perimetry, which has been shown to be a rapid, effective method of detecting glaucomatous visual field loss; (3) Short Wavelength Automated Perimetry (SWAP), which has demonstrated the ability to predict the onset and progression of glaucomatous visual field deficits; (4) The Multifocal Electroretinogram (mfERG) and the Multifocal Visual Evoked Potential (mfVEP), which provide an objective measurement of the visual field. Each of these techniques has presented distinct advantages for the diagnosis and management of glaucoma. *Curr Opin Ophthalmol* 2002, 13:77–84 © 2002 Lippincott Williams & Wilkins, Inc.

Discoveries in Sight Research Labs, Devers Eye Institute, Portland, Oregon, USA.

Correspondence to Chris A. Johnson, PhD, Devers Eye Institute, 1040 NW 22<sup>nd</sup> Avenue, Portland, OR 97210, USA; e-mail: cajohnso@discoveriesinsight.org

**Current Opinion in Ophthalmology** 2002, 13:77–84

## Abbreviations

<b>FDT</b>	Frequency Doubling Technology
<b>mfERG</b>	Multifocal Electroretinogram
<b>mfVEP</b>	Multifocal Visual Evoked Potential
<b>SWAP</b>	Short Wavelength Automated Perimetry
<b>SITA</b>	Swedish Interactive Threshold Algorithm
<b>TOP</b>	Tendency Oriented Perimetry

ISSN 1040–8738 © 2002 Lippincott Williams & Wilkins, Inc.

For many years, quantitative perimetry consisted of determining detection thresholds for a small white target projected onto a uniform white background (white-on-white), using test procedures that were rather time-consuming. The advent of automated static perimetry offered a number of distinct advantages over manual procedures, but the test procedure still consisted of white-on-white threshold determinations performed with time-consuming test procedures. It has only been in recent years that rapid, efficient threshold procedures have become available, along with new perimetric tests that measure visual functions other than white-on-white thresholds. This paper briefly reviews two new efficient threshold test strategies and three new perimetric test procedures, all of which are commercially available.

## Efficient test strategies

### The Swedish Interactive Threshold Algorithms


Swedish Interactive Threshold Algorithms (SITA-Standard and SITA-Fast) are two new test procedures that are available on the Humphrey Field Analyzer II, Model 750 (Humphrey Systems, Dublin, CA). The SITA strategies are based on a forecasting procedure employing Bayesian statistics [1–5]. There are four primary features of the SITA strategies: (1) probability density functions (pdfs), which represent the distribution of thresholds in the population at large for each visual field location (separate pdfs are used to represent normal and glaucomatous populations); (2) likelihood functions (frequency-of-seeing curves); (3) dynamic monitoring of patient response times to interactively pace the test; and (4) comparisons of neighboring locations to adjust the final threshold estimate. Initially, the pdf is used to provide a “best guess” as to the value of threshold. The patient’s response (or lack of response) and the likelihood function are then used to modify the pdf to produce a revised “best guess” as to the value of threshold. This continues until a specified confidence level (or error level in threshold estimate) for the threshold estimate is reached. The difference between SITA-Standard and SITA-Fast is the amount of error that is allowed for the threshold estimate. Details of this type of test procedure may be found in several excellent references [1,3,4,6]. SITA-Standard reduces testing time by 30–50% as compared with conventional Full Threshold procedures, and SITA-Fast produces an additional 30–50% reduction in test time, with the test results comparing very favorably with the findings obtained with Full Threshold procedures [7].

An example of the test results produced by SITA-Standard is presented in Figure 1. The output is quite similar to that obtained for the Full Threshold test procedure, and includes numeric values, gray scale, total and pattern deviation probability plots, reliability indices and visual field indices.

Within the past year, there have been several new findings with regard to the SITA strategies. Sekhar *et al.* [8] confirmed that SITA-Standard and SITA-Fast produced results that were similar to those obtained with Full

Threshold testing, although the average time was approximately 50% less for SITA-Standard and 70% less than SITA-Fast in comparison to Full Threshold testing. Test-retest reliability was excellent for SITA-Standard, but was modestly worse for SITA-Fast. Although the SITA strategies were designed primarily with glaucoma in mind, a recent investigation by Wall *et al.* [9] reports that it is very effective for evaluating other optic neuropathies and hemianopic defects as well. One of the difficulties encountered by clinicians wishing to make the transition from Full Threshold to SITA is how to

**Figure 1. Swedish Interactive Threshold Algorithm-Standard test results**



An example of test results obtained using the SITA-Standard test strategy for the right eye of a glaucoma patient with inferior visual field loss.

compare previous test results obtained using the Full Threshold strategy to new findings obtained using SITA. Although the most desirable alternative is to establish a new baseline with SITA, Heijl *et al.* [10] have derived a method of comparing Full Threshold and SITA results using the Total and Pattern Deviation probability plots when such comparisons are necessary. SITA also recently has been used to perform visual field testing in children [11]. It is clear that SITA has had a significant positive impact on clinical automated perimetric testing.

### **Tendency Oriented Perimetry**

Tendency Oriented Perimetry (TOP) is an efficient test strategy that has been developed for the Octopus 101 and 300 Series perimeters (Interzeag AC, Schlieren, Switzerland) [12]. Tendency Oriented Perimetry uses a staircase procedure, but does so by sequentially evaluating neighboring locations so that each location is tested only once. The visual field is divided into squares of four neighboring test locations and the first location in each square is tested with the initial stimulus intensity. The stimulus intensity for the second presentation is adjusted according to whether or not the patient responded to the first stimulus. The second location in the square is tested with the new stimulus intensity. A similar procedure is employed for the third and fourth members of the square, and a threshold estimate is thus obtained after 4 stimulus presentations. Each threshold estimate is thus based on one direct and three indirect responses from neighbors.

Tendency Oriented Perimetry is approximately four times faster than conventional staircase threshold procedures, needing about three minutes per eye to complete the test [12]. At the present time, there is limited information available concerning the clinical performance of TOP [12]. However, a recent investigation indicates that the short duration of TOP makes it a viable method of performing visual field testing in children down to the age of 6 years [13].

### **Frequency Doubling Technology perimetry**

The frequency doubling effect is created by a low spatial frequency sinusoidal grating ( $<1$  cycle/degree) undergoing high temporal frequency counterphase flicker ( $>15$  Hertz), thereby producing the appearance of twice as many light and dark bars than are physically present (*ie*, the spatial frequency of the bars appears to be doubled). This effect is thought to be mediated by magnocellular (M-cell) retinal ganglion cells with nonlinear response properties. Frequency Doubling Technology (FDT) testing for detection of glaucomatous visual field loss consists of determining contrast sensitivity measures for detecting the frequency doubled stimulus. The commer-

cial version of FDT perimetry (Humphrey Systems, Dublin CA and Welch Allyn, Skaneateles NY, USA) evaluates contrast sensitivity for stimulus detection at 17 or 19 locations throughout the central visual field.


Previous investigations of FDT perimetry indicate that it has high sensitivity and specificity for detection of glaucomatous visual field loss, [14,15] and recent studies have confirmed these findings [16,17,18–20]. In addition, the specificity of the rapid screening test for FDT can be modestly improved by confirming abnormal results [19]. The results of Wu *et al.* [20] suggest that FDT may detect glaucomatous visual field loss earlier than standard automated perimetry. Similarly, Paczka *et al.* [21] found that FDT perimetry had better overall sensitivity and specificity for detection of glaucomatous damage than either scanning laser polarimetry or assessment of retinal nerve fiber layer photographs. Thus, FDT may be sensitive to the earliest changes produced by glaucoma. In addition to its capabilities for detecting glaucomatous visual field loss, FDT perimetry also has been reported to be useful for the detection of certain retinal diseases [22], optic neuritis [23], and other neuro-ophthalmologic deficits [24]. In its present form, however, FDT perimetry has more limited utility for detection of macular disorders [22] and hemianopic defects [24]. Unlike standard automated perimetry, FDT is not affected by the presence of myopia, even in cases of high myopia [25]. Figure 2 presents an example of FDT results obtained in a patient with glaucomatous visual field loss in comparison to findings obtained for standard automated perimetry using the Humphrey Field Analyzer.

### **Short Wavelength Automated Perimetry**

Short Wavelength Automated Perimetry (SWAP) is conducted by projecting a large blue stimulus onto a bright yellow background, which permits the sensitivity of the short wavelength (blue) pathways to be isolated and measured. Short Wavelength Automated Perimetry testing is available on two commercial automated perimeters, the Humphrey Field Analyzer II Model 750 (Humphrey Systems, Dublin CA, USA) and the Octopus 1-2-3 perimeter (Interzeag AC, Schlieren, Switzerland). Previous investigations have established that SWAP is more sensitive than standard automated perimetry for detecting early glaucomatous defects, shows greater progression of existing glaucomatous defects and several other desirable clinical features. An excellent review of the principles underlying SWAP and the results of previous investigations may be found in Sample's recent article [26].

Recent investigations of SWAP have characterized the long term variability of SWAP in comparison to standard automated perimetry, and have established that it exhib-

Figure 2. Frequency Doubling Technology test results



An example of a superior arcuate deficit for Frequency Doubling Technology (FDT) perimetry in comparison to the results obtained for standard automated perimetry with the Humphrey Field Analyzer. Note that the FDT deficit is more extensive than those obtained with standard automated perimetry.

its greater variability [27]. Optimum criteria for detection of glaucomatous visual field loss with SWAP have been developed, and consist of a cluster of four or more points worse than the  $p < .05$  level or a cluster of three or more points worse than the 1% level [28]. An investigation of


methods to detect SWAP progression indicates that using a ranked threshold distribution analysis is closely related to results obtained using pointwise analysis, and is more useful for identifying a diffuse component of change in the visual field [29]. An investigation of a large

group of ocular hypertensive patients revealed that there was a higher prevalence of SWAP deficits than for standard automated perimetry, but that the incidence of new SWAP and standard automated perimetry deficits was equivalent [30]. This indicates that both test procedures are monitoring the same underlying glaucomatous disease process, but that SWAP deficits are detecting the damage at an earlier stage. It is interesting to note that a recent study of SWAP and FDT testing in ocular hypertensives showed that both techniques found deficits that were not picked up with standard automated perimetry,

and that SWAP and FDT abnormalities were highly correlated [31].


Short Wavelength Automated Perimetry deficits have been found to correlate well with glaucomatous neuroretinal rim abnormalities [32], retinal nerve fiber layer defects [33] and optic disc damage [34]. In addition, SWAP has been reported to be superior to standard automated perimetry for identifying patients with progressive optic disc cupping, and may therefore be quite useful for determining progressive changes in glaucoma [35].

**Figure 3. Short Wavelength Automated Perimetry test results**



An example of a superior nasal step that is present for Short Wavelength Automated Perimetry (SWAP) in a patient with evidence of glaucomatous damage to the inferior portion of the optic nerve, but has normal visual fields for standard automated perimetry.



**Figure 4. Multifocal Visual Evoked Potential test results**

Example of correspondence of visually evoked potential perimetry with standard automated perimetry in a glaucoma patient with a superior paracentral arcuate defect. (Reprinted with permission from Klistorner and Graham, *Ophthalmology* 2000, 107:2283–2299).

Figure 3 presents an example of a superior nasal step that was detected by SWAP but is not detected by standard automated perimetry. The patient had evidence of glaucomatous optic disc changes that were consistent with the presence and location of the SWAP deficit, but standard automated perimetry results were within normal limits.

Finally, it has previously been reported that SWAP is also useful in evaluation of ocular diseases other than glaucoma [26]. Recent studies have indicated that SWAP is useful in evaluating diabetic retinopathy [36] and macular function [37].

### Multifocal Visual Evoked Potential

Perimetric assessment of the visual field has traditionally been performed using techniques that require a subjective response on the part of the patient. However, with the development of multifocal electrophysiologic tech-

niques by Sutter *et al.* [38,39], it is now possible to obtain objective measures of visual function from a number of localized regions of the visual field. A detailed description of the principles underlying multifocal techniques is beyond the scope of this review, although there are several excellent references available for readers interested in this information [38–40].

There are two multifocal techniques that are presently being used, the multifocal electroretinogram (mfERG), which measures the local electrical responses of the retina throughout the central (26 degrees radius) visual field and the multifocal visual evoked potential (mfVEP), which measures the localized electrical responses from the primary visual portion of the brain (occipital cortex) for the central (26 degree radius) visual field [40]. Recent investigations of mfERG in patients with glaucomatous visual field loss indicate that they show abnormalities in the mfERG [41–43]. In some in-

stances, glaucoma patients will exhibit reductions in the amplitude of the mfERG, particularly in the most central portion of the visual field [42]. In addition, reductions in a portion of the mfERG waveform, known as the Optic Nerve Head Component, will be reduced or absent in patients with glaucomatous visual field loss [44]. However, these abnormalities in the mfERG can be rather subtle, even in the presence of significant visual field loss for standard automated perimetry. Moreover, the mfERG findings in most glaucoma patients do not correlate well topographically with visual field losses that are present on standard automated perimetry [42]. The usefulness of the mfERG in glaucoma is therefore somewhat limited.

On the other hand, recent investigations of the mfVEP have indicated that it can provide an objective determination of glaucomatous visual field loss that shows very good correspondence with standard automated perimetry results [40,45,46]. Recent investigations within the past year have reported several improvements to the mfVEP technique. Comparisons of mfVEP results with Humphrey Field Analyzer is complicated by the fact that the two tests have different spatial arrangements for their stimulus presentations. In this view, Hood and Zhang [40] developed an interpolation procedure to permit more direct comparisons of mfVEP and Humphrey Field Analyzer test results, and found very good correspondence between the two techniques in terms of characterizing glaucomatous visual field loss. Klistorner and Graham [46] found that by summing responses according to visual field sectors was useful for reducing variability and enhancing waveforms, provided that the sector grouping was conducted appropriately. Klistorner and Graham also found that more robust signals could be obtained by obtaining recordings from four electrode channels rather than one, particularly along the horizontal meridian, and that reproducibility was very good [45].

Figure 4 presents an example of the left eye of a patient with a superior arcuate nerve fiber bundle deficit as measured with standard automated perimetry (top) and the mfVEP (bottom). Note that the mfVEP procedure not only detects the glaucomatous damage, but the topographic location of the mfVEP abnormalities correlate very well with the region of visual field damage for standard automated perimetry. Multifocal visual evoked potential seems to be a promising new technique for objectively assessing the central visual field and will probably be most useful for patients who are difficult to test with conventional perimetric procedures.

## Acknowledgments

The author is a consultant for and receives research support from Welch Allyn (Skaneateles, NY) and Humphrey Systems (Dublin, CA). Supported in part by National Eye Institute Research Grant #EY-03424.

## References and recommended reading

Papers of particular interest, published within the annual period of review, have been highlighted as:

- Of special interest
  - Of outstanding interest
- 1 Bengtsson B, Olsson J, Heijl A, et al.: A new generation of algorithms for computerized threshold perimetry, SITA. *Acta Ophthalmol Scand* 1997, 75:368–375.
  - 2 Bengtsson B, Heijl A.: SITA Fast, a new rapid perimetric threshold test. Description of methods and evaluation in patients with manifest and suspect glaucoma. *Acta Ophthalmol Scand* 1998, 76:431–437.
  - 3 Bengtsson B, Heijl A.: Evaluation of a new perimetric threshold strategy, SITA, in patients with manifest and suspect glaucoma. *Acta Ophthalmol Scand* 1998, 76(3):268–72.
  - 4 Bengtsson B, Heijl A, Olsson J: Evaluation of a new threshold visual field strategy, SITA, in normal subjects. Swedish Interactive Thresholding Algorithm. *Acta Ophthalmol Scand* 1998, 76:165–169.
  - 5 Bengtsson B, Heijl A: Inter-subject variability and normal limits of the SITA Standard, SITA Fast, and the Humphrey Full Threshold computerized perimetry strategies, SITA STATPAC. *Acta Ophthalmol Scand* 1999, 77:125–129.
  - 6 Vingrys AJ, Pianta MJ: A new look at threshold estimation algorithms for automated static perimetry. *Optometry and Visual Science* 1999, 76:588–595.
  - 7 Wild JM, Pacey IE, O'Neill EC, et al.: The SITA perimetric threshold algorithms in glaucoma. *Invest Ophthalmol Vis Sci* 1999, 40:1998–2009.
  - 8 Sekhar GC, Naduvilath TJ, Lakkai M, et al.: Sensitivity of Swedish interactive threshold algorithm compared with standard full threshold algorithm in Humphrey visual field testing. *Ophthalmology* 2000, 107:1303–1308.
  - 9 Wall M, Punke SG, Stickney TL, et al.: SITA Standard in optic neuropathies and hemianopias: a comparison with full threshold testing. *Invest Ophthalmol Vis Sci* 2001, 42:528–537.
  - This investigation compared the results obtained with SITA-Standard and Full Threshold testing on the Humphrey Field Analyzer in patients with non-glaucomatous optic neuropathies and hemianopic deficits. Although the absolute thresholds were slightly higher with SITA-Standard (about 1 dB), the probability plots for the two procedures were highly comparable. The findings indicated that SITA-Standard were at least as good as Full Threshold procedures for characterization of neuro-ophthalmologic visual field loss.
  - 10 Heijl A, Bengtsson B, Patella VM: Glaucoma follow-up when converting from long to short perimetric threshold tests. *Arch Ophthalmol* 2000, 118:489–493.
  - One of the difficulties in switching from one visual field test strategy to another is being able to compare previous visual field results to current findings. In this paper, the authors present a method based on comparison of probability plots in order to compare previous Full Threshold results to SITA for the Humphrey Field Analyzer. An analysis of longitudinal data from 31 glaucoma patients indicates that this approach is feasible and valid.
  - 11 Donahue SP, Porter A: SITA visual field testing in children. *Journal of the American Association for Pediatric Ophthalmology and Strabismus* 2001, 5:114–117.
  - 12 Morales J, Weitzman ML, Gonzalez de la Rosa M: Comparison between Tendency-Oriented Perimetry (TOP) and octopus threshold perimetry. *Ophthalmology* 2000, 107:134–142.
  - This paper introduces a new threshold testing strategy for the Octopus perimetry known as Tendency Oriented Perimetry (TOP). Initial results indicate that the results obtained with TOP compare favorably with findings obtained using conventional staircase procedures.
  - 13 Morales J, Brown SM: The feasibility of short automated static perimetry in children. *Ophthalmology* 2001, 108:157–162.
  - 14 Johnson CA, Samuels SJ: Screening for glaucomatous visual field loss with frequency-doubling perimetry. *Invest Ophthalmol Vis Sci* 1997, 38:413–425.
  - 15 Cello KE, Nelson-Quigg JM, Johnson CA: Frequency doubling technology perimetry for detection of glaucomatous visual field loss. *Am J Ophthalmol* 2000, 129:314–322.
  - 16 Casson R, James B, Rubinstein A, et al.: Clinical comparison of frequency doubling technology perimetry and Humphrey perimetry. *British Journal of Ophthalmology* 2000, 85:360–362.
  - 17 Hollo G, Szabo A, Vargha P: Scanning laser polarimetry versus frequency doubling perimetry and conventional threshold perimetry: Changes during a 12-month followup in preperimetric glaucoma. A pilot study. *Acta Ophthalmologica Scandinavica* 2001, 79:403–407.

- 18 Serguhn S, Spiegel D: Comparison of frequency doubling perimetry and standard achromatic computerized perimetry in patients with glaucoma. *Graefes Archive Clin Exp Ophthalmology* 2001, 239:351–355.
  - 19 Khong JJ, Dimitrov PN, Rait J, et al.: Can the specificity of the FDT for glaucoma be improved by confirming abnormal results? *Journal of Glaucoma* 2001, 10:199–202.
  - 20 Wu L, Suzuki Y, Kunitatsu S, et al.: Frequency Doubling Technology and confocal scanning ophthalmoscopic optic disc analysis in open-angle glaucoma with hemifield defects. *Journal of Glaucoma* 2001, 10:256–260.
  - 21 Paczka JA, Friedman DS, Quigley HA, et al.: Diagnostic capabilities of frequency-doubling technology, scanning laser polarimetry, and nerve fiber layer photographs to distinguish glaucomatous damage. *American Journal of Ophthalmology* 2001, 131:188–197.
- In this study, the authors compared the ability of three diagnostic tests to adequately screen for glaucoma: frequency doubling technology (FDT) perimetry, scanning laser polarimetry (GDx) and nerve fiber layer (NFL) photographs. The findings indicated that FDT had the best diagnostic performance of the three methods evaluated.
- 22 Sheu SJ, Chen YY, Lin HC, et al.: Frequency doubling technology perimetry in retinal diseases: preliminary report. *Kaohsiung Journal of Medical Science* 2001, 17:25–28.
  - 23 Fujimoto N, Adachi-Usami E: Frequency doubling perimetry in resolved optic neuritis. *Invest Ophthalmol Vis Sci* 2000, 41:2558–2560.
  - 24 Thomas D, Thomas R, Muliyl JP, et al.: Role of frequency doubling perimetry in detecting neuro-ophthalmic visual field defects. *American Journal of Ophthalmology* 2001, 131:734–741.
  - 25 Ito A, Kawabata H, Fujimoto N, Adachi-Usami E: Effect of myopia on frequency-doubling perimetry. *Investigative Ophthalmology & Visual Science* 2001, 42:1107–1110.
  - 26 Sample PA: Short-wavelength automated perimetry: it's role in the clinic and for understanding ganglion cell function. *Prog Retin Eye Res* 2000, 19:369–383.
- This is an outstanding review of the underlying basis of short wavelength automated perimetry (SWAP), its clinical utility for evaluation of ocular hypertensive and glaucoma patients, its role in evaluating other ocular disorders and future directions for new research on SWAP.
- 27 Hutchings N, Hosking SL, Wild JM, et al.: Long-term fluctuation in short-wavelength automated perimetry in glaucoma suspects and glaucoma patients. *Investigative Ophthalmology & Visual Science* 2001, 42:2332–2337.
  - 28 Polo V, Larrosa JM, Pinilla I, et al.: Optimum criteria for short-wavelength automated perimetry. *Ophthalmology* 2001, 108:285–289.
- Short Wavelength Automated Perimetry (SWAP) was performed in 95 glaucomatous eyes and 128 normal eyes. The authors then applied various analysis criteria to determine which criteria provided the best sensitivity and specificity for SWAP defects. The findings indicated that a cluster of 4 or more points beyond the 5% probability level or 3 or more points beyond the 1% probability level provided the best performance.
- 29 Kono Y, Sample PA, Emdadi A, et al.: Comparative study between pointwise and ranked threshold distribution analyses of change in serial fields for short-wavelength automated perimetry. *J Glaucoma* 2000, 9:419–427.
  - 30 Demirel S, Johnson CA: Incidence and prevalence of short wavelength automated perimetry deficits in ocular hypertensive patients. *Am J Ophthalmol* 2001, 131:709–715.
- Evaluation of Short Wavelength Automated Perimetry (SWAP) in a large group of ocular hypertensive patients indicated that although the prevalence of SWAP deficits was much higher than for standard automated perimetry, the incidence of new deficits was the same. This suggests that both standard automated perimetry and SWAP are detecting the same underlying glaucomatous disease process, but that SWAP is detecting it earlier.
- 31 Landers J, Goldberg I, Graham S: A comparison of short wavelength automated perimetry with frequency doubling perimetry for the early detection of visual field loss in ocular hypertension. *Clin Experiment Ophthalmol* 2000, 28:248–252.
- A comparison of standard automated perimetry, Short Wavelength Automated Perimetry (SWAP) and Frequency Doubling Technology (FDT) perimetry was performed in 62 glaucoma patients. The results indicated that both SWAP and FDT were able to detect early deficits and may be predictive of future losses for standard automated perimetry. Concordance between SWAP and FDT results was very good.
- 32 Larrosa JM, Polo V, Pablo L, et al.: Short-wavelength automated perimetry and neuroretinal rim area. *Eur J Ophthalmol* 2000, 10:116–120.
  - 33 Mok KH, Lee VWL: Nerve fiber analyzer and short-wavelength automated perimetry in glaucoma suspects: a pilot study. *Ophthalmology* 2000, 107:2101–2104.
  - 34 Ugurlu S, Hoffman D, Garway-Health DF, et al.: Relationship between structural abnormalities and short-wavelength perimetric defects in eyes at risk of glaucoma. *Am J Ophthalmol* 2000, 129:592–598.
  - 35 Girkin CA, Emdadi A, Sample PA, et al.: Short-wavelength automated perimetry and standard perimetry in the detection of progressive optic disc cupping. *Arch Ophthalmol* 2000, 118:1231–1236.
  - 36 Nomura R, Terasaki H, Hirose H, et al.: Blue-on-yellow perimetry to evaluate S cone sensitivity in diabetics. *Ophthalmic Res* 2000, 32:69–72.
  - 37 Remky A, Arend O, Hendricks S: Short-wavelength automated perimetry and capillary density in early diabetic maculopathy. *Invest Ophthalmol Vis Sci* 2000, 41:274–281.
  - 38 Sutter EE, Tran D: The field topography of ERG components in man—I. The photopic luminance response. *Vision Res* 1992, 32:433–446.
  - 39 Baseler HA, Sutter EE, Klein SA, et al.: The topography of visual evoked response properties across the visual field. *Electroencephalography and clinical neurophysiology* 1994, 90:65–81.
  - 40 Hood DC, Zhang X: Multifocal ERG and VEP responses and visual fields: comparing disease-related changes. *Doc Ophthalmol* 2000, 100:115–137.
- This article provides an overview of the use of multifocal ERG and VEP in the evaluation of glaucoma, retinal disease and optic nerve disorders. It also provides guidelines for making appropriate comparisons between multifocal ERG and VEP results and standard automated perimetry findings using the Humphrey Field Analyzer.
- 41 Fortune B, Cioffi GA, Johnson CA, et al.: The relationship between multifocal electroretinogram and standard automated perimetry findings in normal tension glaucoma. In: Weinreb RN, Krieglstein GK, Kitazawa Y, editors. *Glaucoma in the 21st Century*. London: Harcourt, 2000:73–78.
  - 42 Fortune B, Johnson CA, Cioffi GA: The topographic relationship between multifocal electroretinographic and behavioral perimetric measures of function in glaucoma. *Optom Vis Sci* 2001, 78:206–214.
  - 43 Bearse MA, Stamper RL, Sutter EE: Detection of functional abnormalities in glaucoma using a new multifocal ERG paradigm. *Invest Ophthalmol Vis Sci* 2000, 41:S103.
  - 44 Bearse MA, Sutter EE, Sim D, et al.: Glaucomatous dysfunction revealed in higher-order components of the electroretinogram. Washington, D.C.: Optical Society of America, 1996.
  - 45 Klistorner A, Graham SL: Objective perimetry in glaucoma. *Ophthalmology* 2000, 107:2283–2299.
- In this paper, the authors present the results of a comparison study of multifocal VEP and standard automated perimetry in normal observers, glaucoma patients and glaucoma suspects. The findings revealed high sensitivity and specificity for the multifocal VEP as compared to Humphrey Field Analyzer results, as well as good topographic correlations of the location of deficits revealed by the two techniques.
- 46 Klistorner AI, Graham SL: Multifocal pattern VEP perimetry: analysis of sectoral waveforms. *Doc Ophthalmol* 1999, 98:183–196.